Spinbag


IT	Foglietto illustrativo	3
EN	Package leaflet	5
FR	Notice d'utilisation	7
DE	Gebrauchsanweisung	9

4008071900SP

CE 0476

Rev. 01 del 13/10/2023

Simboli / Symbols / Symboles / Symbole

IT Conservare a una temperatura compresa tra +4°C e +30°C
 EN Store between +4°C and +30°C
 FR Conserver à une température comprise entre +4°C et +30°C
 DE Lagerung zwischen +4°C und +30°C

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INTENDED USE

The Spinbag product is an acid concentrate for bicarbonate based hemodialysis and is intended to be used together with a dry sodium bicarbonate cartridge for on-line preparation of hemodialysis, hemodiafiltration or hemofiltration fluids on compatible dialysis systems.

CAUTIONS

The information in this package leaflet must be carefully read and understood before using the Spinbag product.

When the container is empty it shall be discarded according to local regulations.

This package leaflet must always be used together with the operator's manual of the corresponding

dialysis machine which describes the handling instructions.

The operation of the dialysis machine is restricted to persons who are trained in dialysis and who are completely familiar with the operating instructions stated in the operator's manual of the dialysis machine and in this package leaflet.

The concentrate has a low pH. Be careful to avoid concentrate from the Spinbag product splashing into your eyes. However, if you come in contact with the concentrate specific first aid recommendations should be followed, described under FIRST AID.

To be used under medical supervision.

WARNINGS

DO NOT use if any damage is found or suspected.

DO NOT use if the concentrate is not clear and free from visible particles.

DO NOT use after the expiry date written on the bag and on the box labelling.

NOT for infusion.

DO NOT use together with any other type of sodium bicarbonate (NaHCO_3) concentrate than pure sodium bicarbonate concentrate.

DO NOT initiate dialysis without

having followed the instructions in the operator's manual of the corresponding dialysis machine regarding any recommended tests to verify the proper composition of the dialysis fluid.

Use within 24 hours after bag perforation.

Make sure that the composition corresponds to physician's prescription.

Only use on dialysis machines compatible with the Spinbag product and with settings for the actual composition and dilution ratio.

Contact Spindial in case of use with non-compatible dialysis machines, because use of a connection set is required to ensure compatibility.

DO NOT perform any unauthorized modifications or alterations to the product or to the dialysis machine as this may result in a malfunction or have other serious consequences on the safe operation of the dialysis treatment.

FIRST AID**Skin contact**

Remove contaminated clothing. Rinse skin with plenty of water.

Eye contact

Rinse immediately with water for several minutes. Keep eyes wide open. Remove contact lenses, if present and easy to do. Continue rinsing. If irritation remain contact physician.

Ingestion

Wash mouth with water and drink a couple of glasses of water or milk. Do not induce vomiting. If irritation remain contact physician.

DESCRIPTION

The Spinbag product is a polyolefin container containing a concentrated electrolyte solution (e.g. KCl , CaCl_2 , MgCl_2) and acetic acid needed in a bicarbonate based dialysis fluid.

The Spinbag product is available only with glucose.

When the product is attached to the acidic concentrate connector, concentrate is drawn out by the dialysis machine. The concentrate is proportioned by the dialysis machine.

The dialysis machine mixes the acidic concentrate with the sodium bicarbonate concentrate to produce a bicarbonate based dialysis fluid.

Free from bacterial endotoxins.

PVC-free.

SPECIFICATION**Contents**

The Spinbag product is available in two sizes: 3500 ml and 5000 ml. The dilution ratio is 1+44.

Please refer to the table at the end of this leaflet for the exact composition.

Capacity

Bag volume, (ml)	Dilution	Volume dialysis fluid after dilution, (l)
3500	1+44	157
5000	1+44	225

Requirements for mixing

The acetic acid based dialysis concentrate must be mixed with a pure sodium bicarbonate concentrate on a specially adapted dialysis machine to obtain the final dialysis fluid.

Requirements for water

Dialysis machine inlet water quality must comply with local regulations and if no such regulations are available follow ISO 13959.

INDICATIONS

To be used in a bicarbonate based dialysis treatment for patients suffering from acute renal failure, chronic renal failure or acute intoxication with dialysable substances.

CONTRAINDICATIONS

There are no absolute contraindications for bicarbonate based dialysis. The Spinbag product must be used only on the direction of a physician who has evaluated all the pertinent features of the product in relation to the individual patient.

ADVERSE REACTIONS

Some undesirable effects related to dialysis can occur, for instance hypotension, nausea, vomiting and cramps. It is necessary to report any serious incident occurring in connection with the Spinbag device to Spindial and to the competent

authority of the Member State in which the user is established.

INTERACTIONS

The blood concentrations of other drugs can be reduced by the treatment. The blood concentration of certain drugs may need to be monitored and appropriate therapy implemented to correct for removal during treatment. In patients with cardiovascular disease, especially those using cardiac glycoside medications, plasma levels of calcium, potassium and magnesium must be carefully monitored.

INSTRUCTIONS FOR USE

- Remove the protective cap from the valve connector.
- Hang container from eyelet support.
- Fully insert the dedicated connection perforator present on the haemodialysis machine into the access point of the bag.

STORAGE AND SHELF LIFE

Storage

Store between +4°C and +30°C.

Shelf life

Refer to the expiry date stated on the bag and box labelling.

WARRANTY AND LIMITATIONS OF LIABILITY

a) The manufacturer warrants that the Spinbag product has been manufactured in accordance with its specifications and in compliance with Good Manufacturing Practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement only, remedy manufacturing defects in the product becoming apparent within six (6) month from the date of delivery to the customer.

b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability, fitness for a particular purpose, or other warranties, which extended beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the product and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury, or expense arising directly or indirectly from the use of

the product, whether as a result of any defect therein or otherwise.

c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the product, failure or omission to inspect the product before use in order to insure that the product is in proper condition, or any warranty given by independent distributors or dealers.

d) The manufacturer is:
SPINDIAL S.p.A.
Strada Nazionale, 18
43044 Lemignano Collecchio (PR)
Italy

3.5 Litres formulas, dilution 1+44 – concentration expressed in mmol/L

Formula	Na ⁺	K ⁺	Ca ⁺⁺	Mg ⁺⁺	Cl ⁻	CH ₃ COO ⁻	Glucosio
ACS2025	103	2.0	1.25	0.50	108.50	3.0	5.55
ACS2050	103	2.0	1.50	0.50	109.00	3.0	5.55
ACS2075	103	2.0	1.75	0.50	109.50	3.0	5.55
ACS2550	103	2.5	1.50	0.50	109.50	3.0	5.55
ACS3025	103	3.0	1.25	0.50	109.50	3.0	5.55
ACS3050	103	3.0	1.50	0.50	110.00	3.0	5.55
ACS3075	103	3.0	1.75	0.50	110.50	3.0	5.55

5.0 Litres formulas, dilution 1+44 – concentration expressed in mmol/L

Formula	Na ⁺	K ⁺	Ca ⁺⁺	Mg ⁺⁺	Cl ⁻	CH ₃ COO ⁻	Glucosio
ACL2025	103	2.0	1.25	0.50	108.50	3.0	5.55
ACL2050	103	2.0	1.50	0.50	109.00	3.0	5.55
ACL3050	103	3.0	1.50	0.50	110.00	3.0	5.55
ACL3075	103	3.0	1.75	0.50	110.50	3.0	5.55

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